


MRG-059-ZA - Irradiation of [Ru(bpy)3]Cl2 * 6 H2O in the ChemSpeed robot - 2nd pH screening (5.8 -12)

Date: 2024-08-06
Tags: Radiation O2 [Ru(bpy)3]Cl2*6 H2O HTE MRG
Status: Done
Created by: Michael Ringleb

Reaction scheme/sample structure



Literature/reference experiments

Literature	https://doi.org/10.1021/acscatal.6b02595
Reproduction	/
Related experiment	see below AE-334: Preparation of stock solutions for the irradiation of [Ru(bpy)3]Cl2 * 6 H2O AE-JSC-321: Manufacturing and calibration of new 10 mL HTE with sensor spots I

Reagents

Screening of pH

Constant parameters (amounts added for each experiment)

Name	Abbreviation of solution	CAS Number / Experiment Number	Amount [mmol]	concentration used	Molar mass [g/mol]	Volume [mL]	concentration obtained
------	--------------------------	--------------------------------	---------------	--------------------	--------------------	-------------	------------------------

Buffer solution component 1	see varied parameters table below		see varied parameters table below	see varied parameters table below	84.006	see varied parameters table below	/
Buffer solution component 2	see varied parameters table below		see varied parameters table below	see varied parameters table below	105.988	see varied parameters table below	/
Sodium persulfate solution	Ox	7775-27-1	0.051	60 mM	238.09	0.85	6 mM
[Ru(bpy) ₃]Cl ₂ * 6 H ₂ O solution	Ru	50525-27-4	0.085 µmol	40 µmol/L	748.62	0.213	10 µM
Milli-Q water		7732-18-5	/	/	18.015	see varied parameters table below	/

Varied parameters

*final volume of all reaction solutions is 8.5 mL

Experiment	pH value		Abbreviation	end concentration [M]	Volume [mL]	volume Milli-Q water [mL]
MRG-059-ZA-1-1	5.8	Acetic acid/sodium acetate buffer	OAc	0.00835 + 0.0917	0.850	6.375
MRG-059-ZA-1-2	5.8	Acetic acid/sodium acetate buffer	OAc	0.00835 + 0.0917	0.850	6.375
MRG-059-ZA-1-3	5.8	Acetic acid/sodium acetate buffer	OAc	0.00835 + 0.0917	0.850	6.375
MRG-059-ZA-2-1	5.8	kalium hydrogen phosphate/kalium dihydrogen phosphate buffer	KPO4-1	0.0942 + 0.00584	0.850	6.375

MRG-059-ZA-2-2	5.8	kalium hydrogen phosphate/ kalium dihydrogen phosphate buffer	KPO4-1	0.0942 + 0.00584	0.850	6.375
MRG-059-ZA-3-1	8.0	kalium hydrogen phosphate/ kalium dihydrogen phosphate buffer	KPO4-2	0.00652 + 0.0935	0.850	6.375
MRG-059-ZA-3-2	8.0	kalium hydrogen phosphate/ kalium dihydrogen phosphate buffer	KPO4-2	0.00652 + 0.0935	0.850	6.375
MRG-059-ZA-4-1	8.0	ammonia/ ammonium chloride buffer	NH3-1	0.0942 + 0.00575	0.850	6.375
MRG-059-ZA-4-2	8.0	ammonia/ ammonium chloride buffer	NH3-1	0.0942 + 0.00575	0.850	6.375
MRG-059-ZA-5-1	10.0	ammonia/ ammonium chloride buffer	NH3-2	0.0160 + 0.0861	0.850	6.375
MRG-059-ZA-5-2	10.0	ammonia/ ammonium chloride buffer	NH3-2	0.0160 + 0.0861	0.850	6.375
MRG-059-ZA-6-1	10.0	sodium bicarbonate	HCO3	0.070	0.661	6.100
	10.0	sodium hydroxide	OH	0.030	0.464	
MRG-059-ZA-6-2	10.0	sodium bicarbonate	HCO3	0.070	0.661	6.100
	10.0	sodium hydroxide	OH	0.030	0.464	
MRG-059-ZA-7-1	11.0	sodium bicarbonate	HCO3	0.0524	0.4949	5.994
	11.0	sodium hydroxide	OH	0.0476	0.7356	
MRG-059-ZA-7-2	11.0	sodium bicarbonate	HCO3	0.0524	0.4949	5.994
	11.0	sodium hydroxide	OH	0.0476	0.7356	

MRG-059-ZA-8-1	11.0	sodium hydrogen phosphate	NaPO4	0.0859	0.775	6.232
	11.0	sodium hydroxide	OH	0.0141	0.218	
MRG-059-ZA-8-2	11.0	sodium hydrogen phosphate	NaPO4	0.0859	0.775	6.232
	11.0	sodium hydroxide	OH	0.0141	0.218	
MRG-059-ZA-9-1	12.0	sodium hydrogen phosphate	NaPO4	0.0482	0.435	5.989
	12.0	sodium hydroxide	OH	0.0518	0.801	
MRG-059-ZA-9-2	12.0	sodium hydrogen phosphate	NaPO4	0.0482	0.435	5.989
	12.0	sodium hydroxide	OH	0.0518	0.801	

Irradiation Parameters

Power measured using [\[Power Meter\] 843-R-USB + 919P-020-12](#) unless specified otherwise.

Oxygen sensor	Light Source Name	Wavelength [nm]	Power Setting [mW]	Analog input control voltage [V] using Equipment - Joy-it JT-RD6006 DC POWER SUPPLY
FireSting Fiber-Optic Oxygen Meter	Light Source - LCS-0470-50-22	470	75	0.18

Used beam combiner [Name or None]	none
Irradiation distance [cm]	9.5
Thermostat temperature [°C]	/
Stirring speed [rpm]	400
Start time irradiation [s]	see csv/json
End time irradiation [s]	see csv/json

Procedure/observations

Date	Time	Step	Observations
06.08.24		The Experiment AE-JSC-321: Manufacturing and calibration of new 10 mL HTE with sensor spots I was used The protocol for a fully automated workflow (Protocol - Operation of automated workflow for investigation of oxygen evolution) was utilized with a rare earth metal stir bar	used rare earth metal stir bar.jpg ; canula position; utilized vial; position of reagent vials
	17:50-18:10	everything was set up according to the Protocol - Operation of automated workflow for investigation of oxygen evolution	
	18:10-18:17	The python script is loaded as described in Protocol - Operation of automated workflow for investigation of oxygen evolution and the "experiment.yml" is initialized with four experiments (see reagents table for volumes of reactants) with the "run"-parameter = "true" and one at the end of the queue with "run" = "false"	
	18:17	The python script is executed with "python run.py"	
	18:19	AutoSuite program was started	
		from this point on the execution was done fully automatic according to the table with steps executed by the automated platform in the Protocol - Operation of automated workflow for investigation of oxygen evolution	
	18:20	MRG-059-ZA-1-1 was started (OAc placed on place of buffer_solution_1)	septum was pierced 0 times at this point --> rate = see comparison of rates table below
	18:20-18:28	priming of tubings	
07.08.24	6:45		realized that there was a problem in the code --> only one experiment performed before python code did crash because of problem in fitting code --> foto
	6:45-7:30	refilled water bottle, changed part of fitting code (see https://github.com/lamalab-org/photocat-hte/pull/39) --> see foto (code before change; code after changes) emptied vial and rinsed it with fresh water changed lid of vial	

	08:00-08:25	--> Alex prepared new solutions for OAc, Ox and Ru	
	08:30-08:40	--> put new solutions in robot (prepared by Alex Eith (with same concentrations as before))	new solutions could be possible explanation for deviations between MRG-059-ZA-1-1 and 1-2/1-3 but weighing uncertainty below 1%
	08:41	restart of experiments (restart of python code and AuitoSuite app) MRG-059-ZA-1-2 was started (OAc placed on place of buffer_solution_1)	at this point the setpum in the vial lid was pierced 0 times --> rate =see comparison of rates table below
	08:41-08:49	priming of tubings	
	10:05	MRG-059-ZA-1-3 was started (OAc placed on place of buffer_solution_1)	at this point the setpum in the vial lid was pierced 9 times --> rate =see comparison of rates table below
	11:21	MRG-059-Z-2-1 was started (KPO4-1 placed on place of buffer_solution_2)	at this point the setpum in the vial lid was pierced 17 times --> rate =see comparison of rates table below
	12:36	MRG-059-Z-2-2 was started (KPO4-1 placed on place of buffer_solution_2)	at this point the setpum in the vial lid was pierced 27 times --> rate =see comparison of rates table below
	13:45	end of MRG-059-ZA-2-2	at this point the septum was pierced 36 times --> pictures of septum top, bottom
	13:45-13:53	change of lid, replacement of buffer_solution_1 (OAc -> NH3-2), change of experiments in experiment.yml	
	13:53	MRG-059-ZA-5-1 was started (NH3-2 placed on place of buffer_solution_1)	septum was pierced 0 times at this point --> rate =see comparison of rates table below
	14:09	MRG-059-ZA-5-2 was started (NH3-2 placed on place of buffer_solution_1)	at this point the setpum in the vial lid was pierced 9 times --> rate =see comparison of rates table below
	16:25	end of MRG-059-ZA-5-2	at this point the septum was pierced 18 times --> pictures of septum top, bottom
	16:25-16:34	change of lid, replacement of buffer_solution_1 (NH3-2 --> KPO4-2), replacement of buffer_solution_2 (--> NH3-1), change of experiments in experiment.yml	
	16:34	MRG-059-ZA-3-1 was started (KPO4-2 placed on place of buffer_solution_1)	septum was pierced 0 times at this point --> rate =see comparison of rates table below

	17:50	MRG-059-ZA-3-2 was started (KPO4-2 placed on place of buffer_solution_1)	at this point the setpum in the vial lid was pierced 9 times --> rate =see comparison of rates table below
	19:06	MRG-059-Z-4-1 was started (NH3-1 placed on place of buffer_solution_2)	at this point the setpum in the vial lid was pierced 17 times --> rate =see comparison of rates table below
	20:22	MRG-059-Z-4-2 was started (NH3-1 placed on place of buffer_solution_2)	at this point the setpum in the vial lid was pierced 27 times --> rate =see comparison of rates table below
	21:38	end of MRG-059-ZA-4-2	at this point the septum was pierced 36 times --> pictures of septum top, bottom
	21:38-21:47	change of lid, refill of Ox, refill of Milli-Q-water, replacement of buffer_solution_1 (KPO4-2--> HCO3), replacement of buffer_solution_2 (NH3-1 --> OH), change of experiments in experiment.yml	
	21:48	MRG-059-ZA-6-1 was started (HCO3 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	septum was pierced 0 times at this point --> rate =see comparison of rates table below
	23:06	MRG-059-ZA-6-2 was started (HCO3 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	at this point the setpum in the vial lid was pierced 9 times --> rate =see comparison of rates table below
08:08.24	0:24	MRG-059-Z-7-1 was started (HCO3 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	at this point the setpum in the vial lid was pierced 17 times --> rate =see comparison of rates table below
	1:42	MRG-059-Z-7-2 was started (HCO3 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	at this point the setpum in the vial lid was pierced 27 times --> rate =see comparison of rates table below
	3:00	end of MRG-059-ZA-7-2	at this point the septum was pierced 36 times --> pictures of septum top, bottom
	06:50-07:04	change of lid, refill of Ox, replacement of buffer_solution_1 (HCO3 --> NaPO4), change of experiments in experiment.yml	
	07:04	MRG-059-ZA-8-1 was started (NaPO4 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	septum was pierced 0 times at this point --> rate =see comparison of rates table below

	08:22	MRG-059-ZA-8-2 was started (NaPO4 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	at this point the septum in the vial lid was pierced 9 times --> rate =see comparison of rates table below
	09:40	MRG-059-Z-9-1 was started (NaPO4 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	at this point the septum in the vial lid was pierced 17 times --> rate =see comparison of rates table below
	10:58	MRG-059-Z-9-2 was started (NaPO4 placed on place of buffer_solution_1; OH placed on position of buffer_solution_2)	at this point the septum in the vial lid was pierced 27 times --> rate =see comparison of rates table below
	12:16	end of MRG-059-ZA-9-2	at this point the septum was pierced 36 times --> pictures of septum top , bottom

Analysis

Quantitative Interpretation

Comparison of rates

experiment determiner	pH value	utilized buffer	rate constant (utilizing Jacobs fit)	more information
MRG-059-ZA-1-1	5.8	acetate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-1-2	5.8	acetate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-1-3	5.8	acetate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-2-1	5.8	potassium phosphate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)

MRG-059-ZA-2-2	5.8	potassium phosphate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-3-1	8.0	potassium phosphate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-3-2	8.0	potassium phosphate buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-4-1	8.0	ammonia/ammonium chloride buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-4-2	8.0	ammonia/ammonium chloride buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-5-1	10.0	ammonia/ammonium chloride buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-5-2	10.0	ammonia/ammonium chloride buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-6-1	10.0	bicarbonate/hydroxy buffer	0.01845	lid was changed beforehand
MRG-059-ZA-6-2	10.0	bicarbonate/hydroxy buffer	0.01897	septum was pierced 9 times beforehand
MRG-059-ZA-7-1	11.0	bicarbonate/hydroxy buffer	0.02034	septum was pierced 18 times beforehand
MRG-059-ZA-7-2	11.0	bicarbonate/hydroxy buffer	0.02028	septum was pierced 27 times beforehand
MRG-059-ZA-8-1	11.0	sodium phosphate/hydroxy buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-8-2	11.0	sodium phosphate/hydroxy buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)

MRG-059-ZA-9-1	12.0	sodium phosphate/hydroxy buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)
MRG-059-ZA-9-2	12.0	sodium phosphate/hydroxy buffer	NA	for this system no rate was delivered (looks like there was no oxygen evolution at all)

Comparison of [O2]

experiment determiner	pH	utilized buffer	rate constant (from implemented ODE fit from Jacob)	[O2] degassing	[O2] pre-reaction baseline	[O2] reaction	Δ [O2] with underlying oxygen "penetration" through septum	[O2] post reaction baseline	more information
				[μmol/L]	[μmol/L]	[μmol/L]	[μmol/L]	[μmol/L]	
MRG-059-ZA-1-1	5.8	acetate buffer	NA	154 --> 1.5	1.5 --> 39.5	39.5 --> 35 (ca. 1 minute) --> 47 (ca. 9 min)	NA	47 --> 76.6	lid change beforehand
MRG-059-ZA-1-2	5.8	acetate buffer	NA	147 --> 1.3	1.3 --> 7.3	7.3 --> 0.5 (ca. 1 min) --> fluctuation around 0.3 - 0.5 (ca. 9 min)	NA	0.3 --> 7.7	lid change beforehand
MRG-059-ZA-1-3	5.8	acetate buffer	NA	244 --> 1.8	1.8 --> 7.2	7.2 --> 0.5 --> fluctuation around 0.3 - 0.5 (ca. 9 min)	NA	0.3 --> 7.6	septum pierced 9 times beforehand
MRG-059-ZA-2-1	5.8	potassium phosphate buffer	NA	241 --> 1.8	1.8 --> 8.3	8.3 --> 2.5 (ca. 2 min) --> 9 (ca. 8 min)	NA	9 --> 15	septum pierced 18 times beforehand
MRG-059-ZA-2-2	5.8	potassium phosphate buffer	NA	242 --> 1.8	1.8 --> 8.3	8.3 --> 3 (ca. 2 min) --> 9.2	NA	9.2 --> 15.2	septum pierced 27 times beforehand

MRG-059-ZA-3-1	8.0	potassium phosphate buffer	NA	243 --> 1.1	1.1 --> 2.7	2.7 --> 0.1 (ca. 1.5 min) --> 0.9 (ca. 8.5 min)	NA	0.9 --> 1.8	lid change beforehand
MRG-059-ZA-3-2	8.0	potassium phosphate buffer	NA	221 --> 0.8	0.8 --> 2.4	2.4 --> 0.1 (ca. 1.5 min) --> 0.8 (ca. 8.5 min)	NA	0.8 --> 1.6	septum pierced 9 times beforehand
MRG-059-ZA-4-1	8.0	ammonia/ammonium chloride buffer	NA	213 --> 0.5	0.5 --> 0.5	0.5 --> 0 (ca. 10 s) --> 0.3 (ca. 10 min)	NA	0.3 --> 3	septum pierced 18 times beforehand
MRG-059-ZA-4-2	8.0	ammonia/ammonium chloride buffer	NA	214 --> 0.6	0.6 --> 0.5	0.5 --> 0 (ca. 20 s) --> 0.3 (ca. 10 min)	NA	0.3 --> 3	septum pierced 27 times beforehand
MRG-059-ZA-5-1	10.0	ammonia/ammonium chloride buffer	NA	245 --> 0.7	0.7 --> 3.1	3.1 --> 0 (ca. 0.5 min) --> 0.4 (ca. 9.5 min)	NA	0.4 --> 5.2	lid change beforehand
MRG-059-ZA-5-2	10.0	ammonia/ammonium chloride buffer	NA	242 --> 0.7	0.7 --> 3.0	3.0 --> 0 (ca. 0.5 min) --> 0.4 (ca. 9.5 min)	NA	0.4 --> 4.9	septum pierced 9 times beforehand
MRG-059-ZA-6-1	10.0	bicarbonate/hydroxy buffer	0.01845	224 --> 1.4	1.4 --> 22.7	22.7 --> 111 (ca. 2.5 min) --> 129 (ca. 4.25 min) --> fluctuation between 129/130 (ca. 3.25 min)	ca. 90	131 --> 134	lid change beforehand
MRG-059-ZA-6-2	10.0	bicarbonate/hydroxy buffer	0.01897	232 --> 1.8	1.8 --> 3.7	3.7 --> 94 (ca. 3.5 min) --> fluctuation between 93/94 (ca. 1.5 min) --> 82 (ca. 5 min)	ca. 90	82 --> 67	septum pierced 9 times beforehand
MRG-059-ZA-7-1	10.0	bicarbonate/hydroxy buffer	0.02034	218 --> 3	3 --> 4.4	4.4 --> 75 (ca. 2.5 min) --> fluctuation between 74 - 77 (ca. 3 min) --> 68.5 (ca. 4.5 min)	ca. 70	68.5 --> 72.8	septum pierced 18 times beforehand

MRG-059-ZA-7-2	10.0	bicarbonate/hydroxy buffer	0.02028	228 --> 2.9	2.9 --> 4.2	4.2 --> 75 (ca. 2.75 min) --> fluctuation between 74 -76 (ca. 2.5 min) --> 65 (ca. 4.75 min)	ca. 70	65 --> 61.2	septum pierced 27 times beforehand
MRG-059-ZA-8-1	11.0	sodium phosphate/hydroxy buffer	NA	239 --> 1.7	1.7 --> 5.7	5.7 --> 0.8 (ca. 1.25 min) --> 5.8 (ca. 8.75 min)	NA	5.8 --> 9.6	lid change beforehand
MRG-059-ZA-8-2	11.0	sodium phosphate/hydroxy buffer	NA	228 --> 2.4	2.4 --> 6.1	6.1 --> 0.9 (ca. 1.25 min) --> 8.6 (ca. 8.75 min)	NA	8.6 --> 17.4	septum pierced 9 times beforehand
MRG-059-ZA-9-1	11.0	sodium phosphate/hydroxy buffer	NA	225 --> 3.4	3.4 --> 6.5	6.5 --> 0.8 (ca. 1 min) --> 6.7 (ca. 9 min)	NA	6.7 --> 11.1	septum pierced 18 times beforehand
MRG-059-ZA-9-2	11.0	sodium phosphate/hydroxy buffer	NA	222 --> 5.1	5.1 --> 7.3	7.3 --> 0.9 (ca. 1.25 min) --> 9.1 (ca. 8.75 min)	NA	9.1 --> 20.6	septum pierced 27 times beforehand

Conclusions:

- for ZA-1 to ZA-5 and ZA-8,9 no oxygen evolution visible --> potentially no reaction due to low pH or also buffer species
- for ZA-6 to ZA-7 reaction visible --> reaction rates fit to expectations
- duplicates fit well to each other for same experiment (from comparison of [O₂] and reaction rate (were applicable))
- pH 11 seems to be the best pH value for the oxygen production under these circumstances

Linked experiments

- AE-262: Irradiation of PhPDA (AE-257), 1.5 mg/mg SDS, 2 mg/mL PhPDA after 3 d
- AE-265: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- AE-266: Calibration of oxygen and temperature sensor spots in the 10 mL HTE vial
- AE-267: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- AE-271: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, with manual light control
- AE-272: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, with manual light control I
- MRG-059-A: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals
- MRG-059-B: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals II
- MRG-059-C: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals III
- MRG-059-D: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals IV
- MRG-059-E: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals V
- MRG-059-F: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals VI
- MRG-059-Q: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try I
- MRG-059-R: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try II
- MRG-059-S: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try III
- MRG-059-T: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try IV
- MRG-059-U: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try V

- AE-JSC-321: Manufacturing and calibration of new 10 mL HTE with sensor spots I
- AE-323: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- MRG-059-Z: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ in the ChemSpeed robot, fully automated peripherals
- 1st screening of pH (9.2 - 10.4)
- AE-334: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$

Linked items

Equipment - [Irradiation Set-Up](#)

Equipment - [Joy-it JT-RD6006 DC POWER SUPPLY](#)

Light Source - [UHP LED 470 nm](#)

Protocol - [Operation of automated workflow for investigation of oxygen evolution - as of 03.07.2024](#)

Attached files

results_MRG-059-ZA-9-2.png

sha256: 5aa563cf0dd7811e2cadb2890636ee080ed279a7d0a7aee597573e5e3978019d



results_MRG-059-ZA-9-2.json

sha256: 9cc81232d67d628b3f46d2212688aff2f180bc50a412c382e7a0574d8b1554e7

results_MRG-059-ZA-9-2.csv

sha256: 58f86680226c69eba5be100b56b486b1741a10e5cf70e1bf83738d6d16643263

results_MRG-059-ZA-9-1.png

sha256: fc2b6acef2fe6bec1962f3a521a2d68bdd147ca33611e6baa541a2a126395c60



results_MRG-059-ZA-9-1.json

sha256: efcfbdf9e317d02b4631e8a40ef84d1431d315f641dd1486d57330b4ac468c45

results_MRG-059-ZA-9-1.csv

sha256: 90bf7c095912b465981c5b8dbafc78539430eee94b6d9a11532ace8c2eab3081

results_MRG-059-ZA-8-2.png

sha256: 3d78a307d242f258ba7cb82b7cc6a2a191eca3144c3eec184aa833658adc2f4f



results_MRG-059-ZA-8-2.json

sha256: 1f92896751c06444f09f2b917026cdac588cb955c323b330ee864c7d3f7f76f0

results_MRG-059-ZA-8-2.csv

sha256: 2f390dc3967a6334d1c26a5bf994149bfbd83d7f8a40d2d2d25e008a5ff147fb

results_MRG-059-ZA-8-1.png

sha256: 05ae24d40bd7f69167675a3d561864c6cf788fd68aebd8a661dc10c8d27fdc63



results_MRG-059-ZA-8-1.json

sha256: 4713eb85425e2fc43b95fb0234632bcb70ea2d5603093aff174f35b71dd09a50

results_MRG-059-ZA-8-1.csv

sha256: 76f2b8bc12e92cf4ffb3401621542d9c2f99449da39e378f3cde89e67586088a


MRG-059-ZA-8-1and2-9-1and2-lid_top_view.jpg

sha256: 506af4405bb14a75af4de8ef27642f1f92f81baab41bc77d70088099732e2c84



MRG-059-ZA-8-1and2-9-1and2-lid_bottom_view.jpg

sha256: 7d03cff34c465bae5278af2f8e880a9fe8dffa8d3feebc743813e7f5f7f10880



MRG-059-ZA-8-1and2-9-1and2-lid_top_view2.jpg

sha256: 8c7eb9e554001fddcb1bf43656bce00da7280a964a0155147a9451735cc6a3f1



MRG-059-ZA-8-1and2-9-1and2-lid_bottom_view2.jpg

sha256: 9632ae3ddd46ac48de127029d0e4973838f37ba6bd5446d8cf401bad1ef6fd98



results_MRG-059-ZA-9-2.png

sha256: 5aa563cf0dd7811e2cadb2890636ee080ed279a7d0a7aee597573e5e3978019d



results_MRG-059-ZA-9-2.json

sha256: 9cc81232d67d628b3f46d2212688aff2f180bc50a412c382e7a0574d8b1554e7

results_MRG-059-ZA-9-2.csv

sha256: 58f86680226c69eba5be100b56b486b1741a10e5cf70e1bf83738d6d16643263

results_MRG-059-ZA-9-1.png

sha256: fc2b6acef2fe6bec1962f3a521a2d68bdd147ca33611e6baa541a2a126395c60



results_MRG-059-ZA-9-1.json

sha256: efcfbdf9e317d02b4631e8a40ef84d1431d315f641dd1486d57330b4ac468c45

results_MRG-059-ZA-8-2.png

sha256: 3d78a307d242f258ba7cb82b7cc6a2a191eca3144c3eec184aa833658adc2f4f



results_MRG-059-ZA-9-1.csv

sha256: 90bf7c095912b465981c5b8dbafc78539430eee94b6d9a11532ace8c2eab3081

results_MRG-059-ZA-8-2.json

sha256: 1f92896751c06444f09f2b917026cdac588cb955c323b330ee864c7d3f7f76f0

results_MRG-059-ZA-8-2.csv

sha256: 2f390dc3967a6334d1c26a5bf994149bfbd83d7f8a40d2d2d25e008a5ff147fb

results_MRG-059-ZA-8-1.png

sha256: 05ae24d40bd7f69167675a3d561864c6cf788fd68aebd8a661dc10c8d27fdc63



results_MRG-059-ZA-8-1.json

sha256: 4713eb85425e2fc43b95fb0234632bcb70ea2d5603093aff174f35b71dd09a50

results_MRG-059-ZA-8-1.csv

sha256: 76f2b8bc12e92cf4ffb3401621542d9c2f99449da39e378f3cde89e67586088a

results_MRG-059-ZA-7-2.png

sha256: 1e045942577118b454de2251c9b91d9ee52bd59fe2106e47fe2c260a6c7d83ef



results_MRG-059-ZA-7-2.json


sha256: 45c5ed0e3b6dae0f7c2f36dec2e2eb110da9ac6daf096db8e0f0f4b0aebb7b36

results_MRG-059-ZA-7-2.csv

sha256: 3706a222352f7b8a3c016223770068a0c64f688d769d32a6e0a888b41cd3f2b5

fit_MRG-059-ZA-7-2.png

sha256: f9ea785844b6e79302751710b852a1b5037e738fb2c2b24c9ae1141c87649eee



results_MRG-059-ZA-7-1.png

sha256: d34483df4578bc86064c258e5c85d93ab9a8b577452c4c974dd9c0ec748bc2f8



results_MRG-059-ZA-7-1.json


sha256: ff8b0487a1baff11dbb3831b465e06d54e98f15a2b302c67a05423e99f557bcd

results_MRG-059-ZA-7-1.csv

sha256: 1295d5a5289cdec6245496a1b653f970efde6a813eea4109bfe3fa573cf7e229

fit_MRG-059-ZA-7-1.png

sha256: 8c54907b5ceda561f5a1b3198d96625fe372226274171fd717f081d82f4f094e



results_MRG-059-ZA-6-2.png

sha256: df2f252350a48431f785e41ee8c2a492c339e51b53ab1c932eb5a5e8a73af328



results_MRG-059-ZA-6-2.json


sha256: a6e4149acabeedcf7ad38957bffffe8a6d414aef9898dfbe3b57b0afb5413c298

results_MRG-059-ZA-6-2.csv

sha256: 8e75c74112526486f4fd940e5e10ea8499b6e5b337099d335c00536562323711

fit_MRG-059-ZA-6-2.png

sha256: 8b6c192e5a5e26cb3bce2f9fd2c3f820babf8ddf6da1ff02db2a04521702217a



results_MRG-059-ZA-6-1.png

sha256: 8151f58f843b4bd71672977484f6e22793f58bb7efd5c498a47ea3c3b44dede0



results_MRG-059-ZA-6-1.json


sha256: 993de5dabfe9f9946d339a5fabcc502c87ab9c0ccd3e753e5d7fa43dd5b1bbff4

results_MRG-059-ZA-6-1.csv

sha256: 4de0f01a2b852a9ee94e2ea910567c9f7cf70ba078c6b7059f8a2cbf32261a57

fit_MRG-059-ZA-6-1.png

sha256: ee0db064c49a53d86521ed11e78ad3e7ed8db998718f8df8240fcf091e34ad3b



results_MRG-059-ZA-5-2.png

sha256: b550a16f6f529f94ac26dd33a02c73940820f501178099127ec1efda3f237e14



results_MRG-059-ZA-5-2.json

sha256: 1aff90caf49081f8089f9bf8f53d3b0fca89299b0357e13d5a0788fc8ec33f53

results_MRG-059-ZA-5-2.csv
sha256: 9961058b7123a83699d0ae1c084231362874cc3e8426dfb82f20f110f62a2032

results_MRG-059-ZA-5-1.png
sha256: 430ac525fa4458be85286df8d8192502f334c5f3ae4f12743dcad7d0303d07ca



results_MRG-059-ZA-5-1.json
sha256: 3d55d3e8e60751f28f56703d1a300de6b49163ffe8ccbd7de4c3e86d571ade3f

results_MRG-059-ZA-5-1.csv
sha256: d46540555c5ce6ab0312bbb4891fb8b3e96ef86da972b21f7bf8e0f006a4c726

results_MRG-059-ZA-4-2.png
sha256: 4fef3df995b9ac285cb6e04efad8032f53e5f7636b2f0f634ecc62e5c3d10d5f



results_MRG-059-ZA-4-2.json
sha256: c70862aa4345f530558da5d8f125e93cf0e5db016216a3ce5843fcd077c033d1

results_MRG-059-ZA-4-2.csv
sha256: 4ef23a4ec4aa568b8d51841ee68c9171a524ea0f4d5ffa81f464ede9a7feb40b

results_MRG-059-ZA-4-1.png
sha256: 12ca23d8ead21a280b88d5118be0ce0f0abf4d2b852ae846feac48ce3ce46246



results_MRG-059-ZA-4-1.json
sha256: 954a35f4bd4ecbc16a74746fae8b4960ba0f01121d563f4953684591f14294de

results_MRG-059-ZA-3-2.png
sha256: f98306e4833be0128bf662d04c8980ae418122ddb32ee5a9d0c20abee80cc1b2



results_MRG-059-ZA-4-1.csv
sha256: bca1c58126ebec4ae50203e29ed800fcf083853ac095ebed425054cbe143812f

results_MRG-059-ZA-3-2.json
sha256: 6fb55d046e504bbe3ac0fbc7f60f2c5942d303206379cef84468a3cdab915f76

results_MRG-059-ZA-3-2.csv
sha256: b83ca87f21dacfb7a4cb6004e114fd8dc7f42d1a8d1c00511fb630b4f5e391f3

results_MRG-059-ZA-3-1.png
sha256: 4117d8c4560b882423c46a74fd5ec1870ee6721d9a13d6fb5d0e6bfcd7c9cd01



results_MRG-059-ZA-3-1.json
sha256: 7de0154fae9dfab68d54084a38ed295af7e83e2619e3dbf15f8e59409d3a495c

results_MRG-059-ZA-3-1.csv
sha256: d803e4a80ebf512d584dfb3ebf77485058d328408820eaa7fcb8c3841e7f8515

results_MRG-059-ZA-2-2.png
sha256: 99cd8fd482ea6a016c1ca1edc09b7b0d91c45a0797d21c5e96618d4805024713



results_MRG-059-ZA-2-2.json
sha256: c14aaf5bc283e52916e7fcef43a632aa14c7971cbf8caaa0c1e17e3236d22550

results_MRG-059-ZA-2-2.csv
sha256: 2ce42bd660b9e244d5ed367199f9fc7c63c180a201193243abbf4a517419f0f8


results_MRG-059-ZA-2-1.png
sha256: 3f1c0ec760f798e87d288d260429f472fae8dc4a50aecbbfc69ff24ad097a4bf



results_MRG-059-ZA-2-1.json
sha256: d0b5f0a59b07b92b3d245a68e81df6bdff28cbe0a24a727733162a84d4ece992

results_MRG-059-ZA-2-1.csv
sha256: 769bb60ffd767f6ed6276ad2e76895fe25c979ffee31b5bcd1897d45c4a70208

results_MRG-059-ZA-1-3.png
sha256: bcf67317abcdd6e6be4ee2a3c09acf9c2e2bb93d60b61a24f365bf8fc094ac20



results_MRG-059-ZA-1-3.json
sha256: 2fb448cc29098279adf33da38801b226f2ca1f43cac61fcaadd4c3b716b42d7d

results_MRG-059-ZA-1-3.csv

sha256: f6353784a49d5e8b1e465d7398520bb1eceda9c47540a7f714c97ce698b8907f

results_MRG-059-ZA-1-2.png

sha256: 6467569244b5e90997b73b030c499d89970d7be81e6ba7775379dd9d918c7a81



results_MRG-059-ZA-1-1.csv

sha256: 0eb43dcf01dbbe34904746e486c1b2adf7d80856f30e76f4a2fddec5b8598658

MRG-059-ZA-6-1and2-7-1and2-lid_top_view.jpg

sha256: c18dc417fce3122697b3878f216e26e587d0c4cd4368c7199bfd3391948bfff7



MRG-059-ZA-6-1and2-7-1and2-lid_top_view2.jpg

sha256: 77758693e31f4e53954eb78ced1adf3b02816a9ff1a73b52cd805699afc0dfa5



MRG-059-ZA-6-1and2-7-1and2-lid_bottom_view.jpg

sha256: 6cf74c2772b04d76b77261b1a13a4042b7594d8ff7f89f80056cba942b42b9c3



MRG-059-ZA-6-1and2-7-1and2-lid_bottom_view2.jpg

sha256: 04fed10a68242ecae537b11ee4684586dcdeb73a086dc2a944c49948f84a64da




MRG-059-ZA-5-1and2_top_view.jpg

sha256: 3c5cd32051617736216d2d2b63c709a6bbf7b6dfb37fea58636e4487c6636941



MRG-059-ZA-5-1and2_top_view2.jpg

sha256: 8ec77678efdc52f4cac6f7d2e2be29b200ce799f4381a0a0f35328545473587e



MRG-059-ZA-5-1and2_bottom_view.jpg

sha256: d5237feefea3c509599027e6aa1dde40e936212169ac0242c8cdf07d8aa7efb5



MRG-059-ZA-5-1and2_bottom_view2.jpg

sha256: 2196170f794ccdb89985d817ef450845793bc5d21c240197896cc4e79ef38749



MRG-059-ZA-3-1and2-4-1and2-lid_top_view.jpg

sha256: abaa82209493309be3d1a6819b6194e3e2bb0bc068d4a11c1abb6b7f48e93f1b



MRG-059-ZA-3-1and2-4-1and2-lid_top_view2.jpg

sha256: b9a1d46071d85727174d4e3a51da2faaf460ba51f6f56756d359aedadb8d3e37



MRG-059-ZA-3-1and2-4-1and2-lid_bottom_view.jpg

sha256: 5f7da62718ce0f302b57311c13b20880f11d5eacf9ab11add4061fdaf1980777



MRG-059-ZA-3-1and2-4-1and2-lid_bottom_view2.jpg

sha256: f766225ca8d41324654ebddd43233674ba6e41cd26b53908d6fc977f752f2641



MRG-059-ZA-1-2and3-2-1and2-lid_top_view.jpg

sha256: 2aba7b47bd56b00bb790c717bb6ff517c5fe82a5313e7266020e57cc0af344dc



MRG-059-ZA-1-2and3-2-1and2-lid_top_view2.jpg

sha256: 3a18aec491924c8f31c418b22fcca5ae01ad00f0533ff9d0d947df50de6fdd8c



MRG-059-ZA-1-2and3-2-1and2-lid_bottom_view.jpg

sha256: a99ae0f7a611b2b99cd97dd2c1cebd39cd24cb3cdd5bbd0d3bf2fdf1faa5ab84



MRG-059-ZA-1-2and3-2-1and2-lid_bottom_view2.jpg

sha256: 75ff8959cf89e69d8287036b57479a06512e350e2f9f791bb300f16c4a2e9d39



code-before-changes.jpg

sha256: ba0e51133ad656d155d9414adbfc263119c0b454119696b7dbdf78901abf214e



code-after-change.png

sha256: c44e6f8f513faf4ceec62a3e82708276a86362f5d2d8192eb5ffc91ef7721576



error-message-2.jpg

sha256: 86acba0a2ec37c75bd319f66a0ff55ac7c1baec20e68206d831eef38448dd966



Comment

On Jacob Schneidewind wrote:

Notes:

* csv for MRG-059-ZA-1-2 is missing, json for MRG-059-ZA-1-2 and MRG-059-ZA-1-1 is missing



Unique eLabID: 20240805-e9b22ac3b552ec9ea0a1b0d83c22960d7d7b16d5
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1268>