

# MRG-059-G: Accuracy test of liquid handling of SWAVE robot

**Date:** 2024-06-18  
**Tags:** HTE MRG ACCURACY TEST  
**Status:** Done  
**Created by:** Michael Ringleb

## Literature/reference experiments

Literature	<a href="#">Chemspeed manual</a>
Reproduction	/
Related experiment	

## Reagents

Name	CAS Number / Experiment Number	Amount [mmol]	concentration	Equivalents	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Volume [ml]
Deionized water	<a href="#">7732-18-5</a>	/	/	/	/	/	18.015	different volumes (see tables)

## Procedure/observations

Date	Time	Step	Observations
18.06.24		for information regarding parameters of the test, see: "2024_06_18-further_procedure_information.docx"	
	07:30-07:45	all vials which were utilized were weighed on a <a href="#">scale</a> with 0.1 mg accuracy (right scale in ZAF 121) and the masses were entered into a table for storage and later calculation of the deviations. For samples up to 1 mL standard screw top SEC vials were utilized. For samples from 3 mL to 7 mL snap cap glass vials with a volume of 10 mL were utilized.	<a href="#">utilized scale; utilized scale type label</a>

08:24	The rack (80 x 2 mL vial rack) with the vials for the test for Dilutor 1 was placed inside the robot and the program "2024_06_18 - volume_transfer_accuracy_test" was started	overview of setup; overview of program (1); overview of program (2); height of rack overview; height of rak closer look
-08:28	rinsing of needle 1 (Dilutor 1) to flush the tubing -> better accuracy regarding to Chemspeed manual	
-09:45	test for Dilutor 1 (different volumes (see .xlsx for closer information)) vials were always closed with a lid within a timeframe of ca. 5 min after filling during the experiment	
09:45 - 09:50	Rinsing of tubing of Dilutor 2	
09:50 - 11:05	test for Dilutor 2 (different volumes (see .xlsx for closer information)) vials were always closed with a lid within a timeframe of ca. 5 min after filling during the experiment	
11:05 - 11:10	Rinsing of tubing of Dilutor 3	
11:10	start of test for dilutor 3	at low flow rates (1 mL/min --> dilutor sounds trembling)
11:10 - 11:17	vials were closed with a lid directly (0.5 mL) (like for all others before (for Dilutor 1 and 2) --> lunch break -->	
11:17 - 11:42	other vials (0.6 mL, 0.7 mL, 0.85 mL) not directly closed --> stood open until return from lunch break (0.6 mL: ca. 20 min, 0.7 mL: ca. 15 min, 0.85 mL: ca. 5 min)	
11:42	all vials again closed directly after dispensing	
11:58	start of 3 mL aspiration --> problem that the aspiration could not be done because the dilutor could not use the 3 mL/min flow rate for aspiration --> set to 3.5 mL/min --> restarted with rinsing in the beginning (program: "2024_06_18 - volume_transfer_accuracy_test_after_problem")	3.5 mL/min seem to work but nevertheless, the pump/dilutor does not sound great
-13:50	further filling of vials and measuring of them --> see results in "2024_06_18-SWAVE-accuracy handling-measurement values.xlsx"	

## Analysis

Date	Dilutor	Analysis method	Analytical device	Data	Interpretation
18.06.24	1 (1 mL)	gravimetical analysis	<a href="#">analytical scale</a>	in xlsx (sheet with Dilutor 1), also: <a href="#">Dilutor 1_results</a>	for all volumes absolute errors between -2 and -6 $\mu\text{L}$ --> got bigger with increasing volume (however, only 2/3 $\mu\text{L}$ ) relative deviation got smaller with rising volume (from -1.50% (0.2 mL) to -0.5% (1 mL)) absolute deviation increased with increasing volume (from -2 $\mu\text{L}$ (0.2 mL) to -6 $\mu\text{L}$ (1 mL))
18.06.24	2 (1 mL)	gravimetical analysis	<a href="#">analytical scale</a>	in xlsx (sheet with Dilutor 1), also: <a href="#">Dilutor 2_results</a>	or all volumes absolute errors between -4 and -8 $\mu\text{L}$ --> got bigger with increasing volume (however, only 4/5 $\mu\text{L}$ ) relative deviation got smaller with rising volume (from -1.75% (0.2 mL) to -0.75% (1 mL)) absolute deviation increased with increasing volume (from -4 $\mu\text{L}$ (0.2 mL) to -8 $\mu\text{L}$ (1 mL))
18.06.24	3 (10 mL)	gravimetical analysis	<a href="#">analytical scale</a>	in xlsx (sheet with Dilutor 1), also: <a href="#">Dilutor 3_results</a>	larger deviations for small volumes $\leq 1$ mL than with the 1 mL dilutors still errors way below 1% of nominal volume maximum error ca. 40 $\mu\text{L}$ for 7 mL dispensed volume with 7 mL/min flow rate for 5 mL: error ca. 30 $\mu\text{L}$ --> acceptable for reactions (below 1% deviation from desired volume)

## Qualitative Interpretation

### ***For 1 mL syringes:***

generally small deviations (maximum 8  $\mu\text{L}$ ) --> deviations got smaller with increasing volumes (threshold to less than 1% relative deviation = 0.5 mL volume)

deviations for Dilutor 1 smaller than for Dilutor 2 --> use Dilutor 1 in the future for critical volume transfers

in general very good accuracy --> always less than 1% of the nominal volume (1 mL)

### ***For 10 mL syringes:***

deviations bigger than for 1 mL syringes (ca. 1.5 -2 times as high in absolute errors) --> expectable because of bigger nominal volume

relative deviations below 1% for transfers above 3 mL volume

generally always deviation below 1% of nominal volume --> very good

## **Overall:**

all deviations are in negative direction (underfilling of vials) --> probably a systematical error --> however, the error is in such a small range, that we can accept it

accuracies are within the range of Eppendorf-pipettes which would be the manual reference system to use

## **Future recommendations**

### ***Smaller volume transfers (e.g. sacrificial oxidant, buffers, ruthenium solution)***

use for all transfers smaller than 1 mL Dilutor 1 --> best relative and absolute accuracy for this area of transports

use for all transfers smaller than 1 mL 1mL/min as flow rate

let the needle tip dip into the solution to avoid droplet build up and, hence, inaccurate transfers

utilize extra volume and air gap for transfers

### ***Bigger volume transfers (water)***

utilize dilutor 3

utilize slower flow rate (before: 8 mL/min --> now 5 mL/min (for aspiration and dispensing))

include rinsing after water addition

let the needle tip dip into the solution to avoid droplet build up and, hence, inaccurate transfers

utilize extra volume and air gap for transfers

## **Attached files**

2024\_06\_18-further\_procedure\_information.docx

sha256: c475c2d66c3d26cf4f5da28b9d2054ff7eebf0bbcc2d40d04d6f6de2f1f0109a

Chemspeed-Master-Manual.pdf

2024\_06\_18-volume\_transfer\_accuracy\_test\_after\_problem.app

sha256: 60f3e9c238a108c1b1be2e39df6a2694fd124b4956d9498a1f0e2d2b7c5adef6

2024\_06\_18-volume\_transfer\_accuracy\_test.app

sha256: 423cc5b137859875699dd89d94d0fd9715d0d88d38815b581fdeeff1c4cc9467

2024\_06\_18-SWAVE-accuracy-handling-measurement-values.xlsx

sha256: 54ed805c0440b66298a97db91933e32eeb129e6b6a4317d458a6934e09d91d3d


Dilutor-3.png

sha256: 7ce0d6cc4c661006dd2d9420af55674e540c31e5780b70f8ed3f6578d35a214a



Dilutor-2.png

sha256: 1b635855c8d15a7fe74f3dbd0fd6d74b2229c5150a28a8dc367a6727326c9f21



Dilutor-1.png

sha256: 9445f49767b249da9a57445a3804d6f16fa1889489ac3765f303d30e3641a8cb



height\_of\_tray2.jpg

sha256: e5df2a0df11f1f17d4a41c4b390b765404dc1fc9265a98158d4e8a0bd08b847e



height\_of\_tray1.jpg

sha256: 6a4121fb80f564f23f3f90eda422e1a167a9a73acf53840e62c1177af61576f1



program\_for\_accuracy\_test\_overview2.jpg

sha256: baf60164cf4c55ff7bb0ad3a01745537212fbadf9b23435a88351346a1e0ac9c



program\_for\_accuracy\_test\_overview1.jpg

sha256: 9bb338e866820a0cc136a90b96608bb552b407180761c8db18ad3afa13601f2a



overview\_of\_setup.jpg

sha256: 0429bf2397e2ae362644a69d5ebd690bd6cf7071db9e5a8f8ffd0df63c060270



Utilized\_scale\_type\_label.jpg

sha256: e4d87dee09a152993a9dad9796ec666375b254b0183d880a3bfd5ef15e82deba



utilized\_scale.jpg

sha256: 702fd07d3f87d42240d2017eb8bb0adb622be7b7e041ed2d4221dacc6f96a3e1



2024\_06\_18-volume\_transfer\_accuracy\_test\_after\_problem.app

sha256: 60f3e9c238a108c1b1be2e39df6a2694fd124b4956d9498a1f0e2d2b7c5adef6

2024\_06\_18-volume\_transfer\_accuracy\_test.app

sha256: 423cc5b137859875699dd89d94d0fd9715d0d88d38815b581fdeeff1c4cc9467



Unique eLabID: 20240618-44ce5758ddfc18560fd2ef13a012e058501d5a97  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1111>