

AE-413: Titration of Na₂CO₃ with NaHCO₃ and vice versa

Date: 2025-01-08

Tags: pH AE Titration Calibration pH meter

Category: HTE

Status: Done

Created by: Alexander Eith

Reaction scheme/sample structure

Chemical Formula: CNa₂O₃

Molecular Weight: 105,98754

NaHCO₃

Na₂CO₃

Na₂SiF₆

Chemical Formula: CHNaO₃

Molecular Weight: 84,00577

Chemical Formula: F₆Na₂Si

Molecular Weight: 188,05496

Literature/reference experiments

Literature	/
Reproduction	/
Related experiment	AE-345: Titration of sodium hydrogen carbonate with Na₂SiF₆ Organisational - AE-412: Calibration of pH-Meter

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]
NaHCO ₃	144-55-8	/	15	1260.1	1259.23	84,007	/	/
Na ₂ CO ₃	497-19-8	/	15	1589.8	1588.24	105.99	/	/
H ₂ O	/	/	/	/	/	18.02	0.998	300

Procedure/observations

Date	Time	Step	Observations
07.01	14:30	NaHCO ₃ and Na ₂ CO ₃ were weighted into seperate 15 mL snap on cap vials	

	15:15	To both water was added (approx. 15 mL each)	
		The vial were shaken and vortexed (Equipment - VWR® VV3, Vortex Mixer stage 4/6), till everything dissolved	
		The obtained solution were transferred into a 100 mL and 50 mL volumetric flask for each. The solution was splitt approx. 2 to 1 between the 100 mL and 50 mL flask	
		All volumetric flask were filled with Milli-Q water	
	- 15:25	The NaHCO ₃ solutions were combined in a 250 mL beaker, the Na ₂ CO ₃ solutions were combined in a 150 mL beaker	For both 0.1 M solution were obtained
		Generel workflow for the titration: In a 50 mL snap on cap vial 10 mL of one solution was added using a 5 mL volumetric pipett. The solution was stirred at 550 rpm and the pH electrode was submerged. The other solution was added using a 25 mL burette (for exact amounts and recoded pH values see Analysis table)	
	15:40	The titration was done with Na₂CO₃ in the burette and NaHCO₃ in the vial	-1 setup.jpg
	15:55	The titration was done with Na₂CO₃ in the burette and NaHCO₃ in the vial	-2
	16:10	The titration was done with NaHCO₃ in the burette and Na₂CO₃ in the vial	-3
	16:15	The titration was done with NaHCO₃ in the burette and Na₂CO₃ in the vial	-4

Analysis

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Processed Data	Interpretation
------	------	-------------	-----------------	-------------------	---------	----------	----------------	----------------


10.09	11:30	AE-345-1	pH	AE-344: Calibration of pH-Meter	H2O	raw data.jpg	AE-413.xlsx, AE-413-for_python.CSV	See Results or excel file
	13:20	AE-345-2	pH	AE-344: Calibration of pH-Meter	H2O	raw data.jpg	AE-413.xlsx, AE-413-for_python.CSV	
	14:05	AE-345-3	pH	AE-344: Calibration of pH-Meter	H2O	raw data.jpg	AE-413.xlsx, AE-413-for_python.CSV	
	14:40	AE-345-4	pH	AE-344: Calibration of pH-Meter	H2O	raw data.jpg	AE-413.xlsx, AE-413-for_python.CSV	

Results

Looks ok.

Example plot (obtained with AE-413-for_python.CSV): Figure_1.png

Used for plotting and fitting of power function: AE-413-c.py, gives good fit, first 20 data points are omitted



For calcaultion of ratio(NAHCO3/Na2CO3) at given pH:

AE-413-calc.py

can be used in range of 8.6 to 10.9, extreme values (pH = 8.4 and 11.1 can also be used)

Fit in good agreement (max. error of less than 0.1 pH units) compared to <https://www.aatbio.com/resources/buffer-preparations-and-recipes/carbonate-bicarbonate-buffer-ph-9-2-to-1>

Linked experiments

- [AE-344: Calibration of pH-Meter](#)

- [AE-345: Titration of sodium hydrogen carbonate with Na₂SiF₆](#)

Organisational - [AE-412: Calibration of pH-Meter](#)


Linked item

Equipment - [VWR® VV3, Vortex Mixer](#)

Attached files

AE-413-c-plot.png

sha256: 926096cd039129f6377b4b6ee8cd1cff533e9e7bb4a31b238fde383cf714649c



AE-413-calc.py


sha256: 6f8104a8eb43620550c409a4b8567639fac76ba90657aedd737599215125d684

AE-413-c.py

sha256: 01fd24cdb95c11c46e680273b955b51b132dc869fa6e36e81255d98dde862c15

Figure_1.png

sha256: 678f8a178fc22e1edd9ca38d7e52e4593c7cd09aa45e916fc53f944fbe970738



AE-413-a.py

sha256: 653fa8bfdd76efecbc1448ae31947a5e49760ef3111434fb9730492a22b9b824

AE-413-for_python.CSV

sha256: a0d14b79478552e15a6973e009fe5db5e0eb3a44cca47ecd2d687dffc42eaab8

AE-413.xlsx

sha256: 6bfb630150b4d7a84a0ba960a76a5f3fb6009ace31620940e7f92c8a8c052d08

setup.jpg

sha256: 70ff7d959cb49cb2cb9ea6c8db60ab88a835f9173245f90e498f7491745360ce



raw-data.jpg

sha256: 8e98e61482c17fe557adcacd29bc8e57a4af8b25f605de4dc097a12fdb61858



Unique eLabID: 20250108-4279404eae193749d3abb89950f86a8696efe672
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1613>