

Operation of automated workflow for investigation of oxygen evolution - as of 03.07.2024

Date: 2024-06-28
Tags: HTE Irradiation MRG
Category: Protocol
Created by: Michael Ringleb

Goal

Description of the workflow necessary to operate the automated setup for the investigation of oxygen evolution.

Prerequisites and preparation

1. Preparation of the solutions necessary for the planned reaction set. Also note the necessary volumes of the solution for each reaction (necessary for later input into "experiments.yml").

E.g.:

| Name | CAS Number / Experiment Number | Amount [mmol] | concentration | Equivalents | Mass _{theo} [mg] | Mass _{exp} [mg] | Molar mass [g/mol] | Volume [ml] |
|------------------------------|---|---------------|---------------|-------------|---------------------------|--------------------------|--------------------|-------------|
| Sodium bicarbonate | Experiment - AE-270: Preparation of stock solutions for the irradiation of [Ru(bpy)3]Cl2 * 6 H2O (AE-270-1) | 0.868 | 1.09 M | / | / | / | 84.006 | 0.8 |
| Sodium carbonate (anhydrous) | Experiment - AE-270: Preparation of stock solutions for the irradiation of [Ru(bpy)3]Cl2 * 6 H2O (AE-270-2) | 0.632 | 0.79 M | / | / | / | 105.988 | 0.8 |
| Sodium persulfate | Experiment - AE-270: Preparation of stock solutions for the irradiation of [Ru(bpy)3]Cl2 * 6 H2O (AE-270-3) | 48 | 60 mM | / | / | / | 238.09 | 0.8 |
| [Ru(bpy)3]Cl2 * 6 H2O | Experiment - AE-270: Preparation of stock solutions for the irradiation of [Ru(bpy)3]Cl2 * 6 H2O (AE-270-4) | 0.24 µM | 0.3 µM | / | / | / | 748.62 | 0.8 |
| Milli-Q water | 7732-18-5 | / | / | / | / | / | 18.015 | 5.3 |

2. Determine the parameters under which the reaction shall be conducted:

E.g.:

Power measured using [\[Power Meter\] 843-R-USB + 919P-020-12](#) unless specified otherwise.

| Oxygen sensor | Light Source Name | Wavelength [nm] | Power Setting [mW] | Analog input control voltage [V] using Equipment - Joy-it JT-RD6006 DC POWER SUPPLY |
|--|---|-----------------|--------------------|---|
| FireSting Fiber-Optic Oxygen Meter | Light Source - LCS-0470-50-22 | 470 | 75 | 0.18 |

| | |
|-----------------------------------|------|
| Used beam combiner [Name or None] | none |
| Irradiation distance [cm] | 9.5 |
| Thermostat temperature [°C] | / |
| Stirring speed [rpm] | 400 |

3. Think about a determination scheme for the reactions which are performed. (E.g.: First reaction = MRG-055-A, Second = MRG-055-B...)

Steps for user to do

| Step number | Step description |
|-------------|---|
| 0 | Get the firesting, connect it to the Laptop and load the configuration: "240717_Instrument_Settings_HTE_Vial_10mL_Liquid_Phase" |

| | |
|-----|--|
| 0.1 | Connect firesting to group laptop (USB connection) and open PyroWorkbench on laptop (ignore the update prompt) |
| 0.2 | Goto Instrument-->Firesting O2 (2 Channels) (A) --> Load instrument configuration --> FireSting Data --> 240717_Instrument_Settings_HTE_Vial_10mL_Liquid_Phase (until 15.10.24) or 241015_Instruments_Settings_HTE_Vial_Liquid_Phase (from 16.10.24) or 241113_Instruments_Settings_HTE_Vial_Liquid_Phase (from 18.11.24) or xxxx.x.x.x.x.. (from xxx) -->Load |
| 0.3 | Disconnect Firesting from the laptop --> ready to go |
| 1 | Setup has to be installed. See: Overview-of-connections-in-setup.jpg , do not connect the magnetic stirrer's plug to the electrical outlet just yet |
| 2 | Switch the programmable power supply on and also press the "on" button to start it (initial voltage should be 0 V) |
| 3 | Switch the lamp power supply (13 A) on and place the lever on the front to "Trigger mode", make sure once more, that it i plugged into the right power outlet (slave outlet of the master/slave power outlet(white)) |
| 4 | Open the argon valve on the wall and check if the reservoir bottle is filled with Milli-Q-water to the top |
| 6 | Placing necessary solutions in respective Chemspeed 13 mL vials |
| 7 | Placing Chemspeed vials at the respective places of the stock solution rack (see AutoSuite program for details about the position). |
| 8 | Placing a rare earth metal stirring bar in the reaction vial, closing it with a slit septum lid (fully close the lid (screw it in as hard as possible)) |
| 9 | Also bend the degassing needle (small blue canula) and disconnect it from the luer adapter from the argon line. Then insert it to the vial so that the canula tip inside the vial is slightly above the level of the sensor spots (to ensure that no backflow from the vial into the canula happens during filling of the vial) (see photo). Furthermore, the needle should be inserted through the septum not close or in the slit area: Not like this but rather in this way (close to the edge of the septum between the lines of the star-shaped puncture in the septum) (Degassing needle position like this) to avoid high leakage rates evoked by septum collapse. At last, reconnect the luer adapter of the argon line back to the canula. Also check to NOT scratch the sensor spots with the canula. |
| 10 | Then insert the vial to the vial holder (has to clip in) |
| | Connect the plug of the magnetic stirrer to the elctrical outlet and check if the stirring speed on the stirrer is set to 400 rpm |
| 11 | Open python program (click on Visual studio Code symbol in the taskbar) and defining the experiments which should be performed in "experiment.yml" (click on experiment.yml and enter the values for the experiments("buffer_solution_1" = sodium carbonate, "buffer_solution_2" = sodium bicarbonate). After four experiments with parameter run : true, add one with run : false (necessary to change the lid later on) --> save the yml (ctrl+s) |
| 13 | open Anaconda on the IPC of the robot |
| 14 | in Anaconda --> Environments --> base(root) --> press green triangle --> Open Terminal. |

| | |
|--------|--|
| 15 | in Terminal: change directory to folder with run.py (recently: "cd Documents"--> "cd photcat_hte" --> "cd photocat-hte" --> "cd experiments" |
| 16 | in Terminal: execute run.py ("python run.py") |
| 17 | Python and arduino (for degassing valve and lamp outlet) code available at: https://github.com/lamalab-org/photocat-hte |
| 18 | Open Chemspeed AutoSuite program (Desktop\Michael\HTE photocatalysis project\ tests for writing of csv from AutoSuite\ "2024_08_02 - final_version_slower_emptying" (from 02.08.24) or "2024_07_03 - final_version_without degassing valve in autosuite"(until 02.08.24)), click the green triangle in the upper part ("Run Application in Executor in Real) --> the AutoSuite Executor is opened --> click start in the Executor program |
| 19 | Close the lid of the robotic platform but make sure to not squash/squeeze any electrical or inert gas lines! |
| 20 | From this point onwards, the automated setup will perform all steps required for the formulation of a photocatalytic reaction solution and it's activity test (for steps: see table below) |
| 21 | When the first degassing of a cohort (four experiments) is performed: Check if there are bubbles visible which indicate the flow of argon through the solution --> should normally be the case, just check to be extra sure |
| 22 | Turn the light in the room on and place signs which indicate that it should be left on (to ensure same experimental conditions for every conducted experiment) |
| 23 | After four experiments have been performed: Switch the lid of the vial for a new one (make a foto of the old lid and save it to the experiment) --> also screw this new one in firmly |
| 24 | Replace the experiments in "experiment.yml" with a new set of experiments, save them (ctrl + s) and execute "python run.py" again in the Terminal window of the Anaconda terminal |
| 24.1 | Check if the tubing between the 4-needle head and the syringe pumps contains bubbles |
| 24.1.1 | If it contains bubbles: Exit the dialogue window shown in the AutoSuite Executor window by clicking "ok" --> directly afterwards stop the AutoSuite Executor by pressing the "stop" button also stop the python script and restart it (to load the correct experiment) Click on the first macrotask ("Prime tubing ...") in the AutoSuite Executor and then restart the program by clicking onto the "Start" button (the one in the middle with the green triangle) |
| 24.1.2 | If it does not contains bubbles: Click the "ok" button in the dialogue window shown in the AutoSuite Executor --> the program will automatically continue |
| 24.2 | --> do prior points as often as necessary for the planned set of experiments (e.g. if you want to do 12 reactions in total you have to change the values in the "experiment.yml" twice) --> always close the lid of robotic platform again after changing the vial lid. |

| | |
|----|--|
| 25 | After all planned reactions have been performed: Stop python script (if it is still running (ctrl, shift, c at the same time)), Stop the program in the AutoSuite Executor, Disassemble the setup, Clean everything and bring it back to the right place |
| 26 | Close the argon valve on the wall |
| 27 | Export all result logs from the Chemspeed's IPC (path: "Documents\photocat_hte\photocat-hte\experiments\logs") (USB-stick necessary) and save the csv/json files to the elab experiment |
| 28 | Evaluate the reaction outcomes and note these |


Steps in the automated protocol executed by the Chemspeed platform

see pictures of process below ([picture 1](#), [picture 2](#))

| Step number | Step description | approximate time required |
|-------------|---|--|
| 0 | If it is the first reaction of the day , the robot will rinse the required tubings to ensure that they are free of air bubbles and potential impurities. | 7 min |
| | If it is not the first reaction of the day , the robot will directly start with the points below | |
| 1 | Dispensing respective amounts of Milli-Q water, buffer component 1 (e.g., sodium carbonate), buffer component 2 (e.g., sodium bicarbonate), sacrificial oxidant and ruthenium catalyst solution (in this order) into the reaction vial while the emerging reaction solution is stirred at 400 rpm | ca. 25 min |
| 2 | afterwards the solution is automatically deoxygenized by argon sparging | recently 20 min (depending on user input to python code) |
| 3 | then, the solution is allowed to equilibrate in order for the oxygen levels in the headspace and solution to balance | recently 15 min (depending on user input to python code) |
| 4 | afterwards the solution is irradiated with a 470 nm LED to start the reaction | recently 10 min (depending on user input to python code) |
| 5 | after the end of irradiation, again an equilibration time is given | recently 5 min (depending on user input to python code) |

| | | |
|---|--|--|
| 6 | afterwards, the vial is automatically emptied, washed with 9 mL of Milli-Q water and emptied again to prepare it for the next round of experiments | 7 min |
| | | total time for 1 experiment (ca. 1 h 20 min - 1 h 30 min) |

Reaction scheme/sample structure



Literature/reference experiments

| | |
|------------|---|
| Literature | https://doi.org/10.1021/acscatal.6b02595 |
|------------|---|

Comparison to previous protocol

| Category | Old protocol | New protocol |
|--------------|--|--|
| Degassing | different outcome because septum was sealed differing strong | needle without tubing is inserted into the septum to release excess pressure even if septum is very tightly sealed |
| Septum lid | no instruction | tighten as much as possible |
| Degassing II | Green canula utilized --> some problems with bending, etc. | blue canula utilized --> should be pre-bent before inserted to vial (addition of instruction to be careful with touching the sensor spots with the canula) |

| | | |
|--------------------------------|---|--|
| Degassing III | Degassing valve control was previously included in the AutoSuite as RS232 communication device --> had the problem that the port was blocked for further communication (e.g. to switch the master/slave outlet) | Deleting of RS232 device from AutoSuite app --> Inclusion to the python code (sending of "1" (string) for opening of degassing valve and "0" for closing) |
| Stirring | 500 rpm --> problems with recovering of spinning of stirrer after emptying of vial | 400 rpm |
| Guidance for glas fibre cables | no guidance --> much bending pressure on the holders for the glass fibre probes | made holes in box at height of glass fibre probe holders --> no bending pressure to be expected |
| Control of lamp power supply | Problems with small "lightnings" (probably due to short circuit in switch on power supply (Trigger <-> Manual)) | Introduction of master/slave outlet to control the power connection to the lamps power supply --> utilizing a USB elongation cable and a relay to switch the (inspiration source) --> picture of setup Python code sends "ON" or "OFF" to Arduino to switch the relay and hence turn power outlet on or off |
| degassing canula position | degassing canula was forced until the bottom of the vial --> potential problems with scratching of the sensor spots | canula only inserted until ca. 1/3 of the vial (looking from top of lid) --> less risk of scratching |
| | | |

Change from "2024_07_03 - final_version_without degassing valve in autosuite.app" to "2024_08_02 - final_version_slower_emptying.app"

- the emptying of the vials was decreased in velocity
- change of aspiration speed for aspirations from vial after end of post-reaction baseline from 20 mL/min to 10 mL/min
- hope that with that all the liquid gets taken up from the vial, avoiding rests of the previous reaction or rinsing water are carried over to next reaction

Linked experiment

- [AE-JSC-321: Manufacturing and calibration of new 10 mL HTE with sensor spots I](#)

Linked items

Archived - [Protocol - Operation of automated workflow for investigation of oxygen evolution - as of 03.06.2024 until 18.06.2024](#)

Archived - [Protocol - Operation of automated workflow for investigation of oxygen evolution - as of 18.06.2024](#)

Archived - [Protocol - Operation of automated workflow for investigation of oxygen evolution - as of 27.06.2024 \(from MRG-059-M\)](#)

Equipment - [Irradiation Set-Up](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 1\)](#)

Equipment - [Joy-it JT-RD6006 DC POWER SUPPLY](#)

Light Source - [UHP LED 470 nm](#)

Attached files

241212_Instruments_Settings_HTE_Vial_Liquid_Phase.ini
sha256: a6cd38c0d11c2606bb3ab08c0aebb16c367a4ac287c2b8893c5fc9bd84e233c0

needle_depth.jpg
sha256: ba4f6b67750774598747421ff8ab0971bd34d2704ee4c202daebb411731fd04b



2024_08_02-final_version_slower_emptying.app

sha256: 6996a62e7a8f15af2ef8509493edfd32f0c4f0699d771d71a8f1d9f61fe3d0c6

2024_07_04-Degassing-valve-configuration-and-all-connections-scheme.jpg

sha256: 2937c9e482402bda00bacd46dd5bac44c659a5a825867b8fbd68de7f0bcd550b



2024_07_03-final_version_without-degassing-valve-in-autosuite.app

sha256: 53a3f7391efa6d6b743e1c4bfd73dbe90bb1499c4b36d1b672e923fae1f52f7

degassing_canula_depth.jpg

sha256: 1e4c7333ede9f68012338fd1b80a3a16abaa043943b5dd38293d7b33845c8da5



program-overview2.jpg

sha256: 2074552bcfa0c527d3c86e21d79483ede43998aaa3e7ace378fe3ccdb18722dc



program-overview.jpg

sha256: c6b4d7570c3d2b1de741c2cf047c46516c76f8e466612d5161b2aa84c0412e3e



degassing-valve-relay.jpg

sha256: fce4efab9e1b1e5036ef7bde1cb2b991c3aac6f25c28d6a5dffe950bd577b0f2



Power-outlet-relay.jpg

sha256: 6e3e58c8dec70f80356323f579c789ddc2ab1b4346448cc50e126631cc873817




2024_03_07-lamp_outlet_and_degassing_sketch.ino

sha256: cd0166d6a576ba14a0ae188800c34901e1938484102aea5dc2036456cfadeaec

vial-position-in-holder.jpg

sha256: f7c07b72756b50bc2745ffd3225c9e8d1a4b1d32270d0e83ae72f1dc46bbadff



good-needle-position.jpg

sha256: 7fa43daf2d4ee6b5a0d0ff3279b9a1592ce1b08f3a3991fcfbf2dccd1c36c339



canula-to-be-utilized.jpg

sha256: 8f9f7d9a8afc20944e452421be0310d80ef5500e2e0b4bae3e8da537df2c6d60



Comment

On Alexander Eith wrote:

* Adding, that master/slave socket needs to be turned on



Unique eLabID: 20240628-04a1e9200dd1b0ed41ff6a478b60c94d3933629f

Link: <https://elab.water-splitting.org/database.php?mode=view&id=111>