

MRG-059-ZP: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - repetition of -ZN-20 and -ZO-7

Date: 2025-02-03


Tags: Radiation O2 [Ru(bpy)₃]Cl₂*6 H₂O AE reference analytics HTE MRG

Category: HTE

Status: Done

Created by: Michael Ringleb

Reaction scheme/sample structure



Literature/reference experiments

Literature	https://doi.org/10.1021/acscatal.6b02595
Reproduction	/
Related experiment	<p>HTE - AE-434: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O</p> <p>HTE - AE-404: Manufacturing, Test and calibration of HTE vial with O₂ and T spot</p> <p>HTE - AE-406: MRG-059-ZM: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - reproduction of standard conditions - new vial</p> <p>HTE - MRG-059-ZN: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - screening of pH, [Ru], [Ox] and Irradiance</p> <p>HTE - MRG-059-ZO: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - repetition of experiments from -ZN with not reproducible results</p>

Reagents/Parameters table

for content of the stated solutions see [HTE - AE-434: Preparation of stock solutions for the irradiation of \[Ru\(bpy\)₃\]Cl₂ * 6 H₂O](#)

for calculations see:

Run	Used [Ru] solution	V [Ru] [mL]	obtained conc. [μM]	Used [Ox] solution	V [Ox] [mL]	obtained conc. [mM]	Used [NaHCO3] solution	V [NaHCO3] [mL]	obtained conc. [M]	Used [Na2CO3] solution	V [Na2CO3] [mL]	obtained conc. [M]	obtained pH according to HTE - AE-413: Titration of Na2CO3 with NaHCO3 and vice versa	V [H2O] [mL]	power setting [V]	obtained power output [mW] AE-341: Power measurement with different power settings from Joy-it with 470 nm
1	Ru II	0.425	10	Ox III	0.425	6	NaHCO3 II	0.753	0.0886	Na2CO3 II	0.324	0.0114	9.2	6.573	0.18	75
2	Ru II	0.425	10	Ox III	0.425	6	NaHCO3 II	0.602	0.0708	Na2CO3 II	0.828	0.0292	9.2	6.221	0.18	75

*final volume of all reaction solutions is 8.5 mL

Irradiation Parameters

Power measured using [Power Meter] 843-R-USB + 919P-020-12 unless specified otherwise.

Oxygen sensor	Light Source Name	Wavelength [nm]	Power Setting [mW]	Analog input control voltage [V] using Equipment - Joy-it JT-RD6006 DC POWER SUPPLY
FireSting Fiber-Optic Oxygen Meter	Light Source - LCS-0470-50-22	470	see table above	see table above

Used beam combiner [Name or None]	none
Irradiation distance [cm]	9.5
Thermostat temperature [°C]	/
Stirring speed [rpm]	400
Start time irradiation [s]	see csv/json
End time irradiation [s]	see csv/json

Procedure/observations

Experiment was done by AE

Date	Time	Step	Observations
03.02	16:15	The Calibration HTE - AE-404: Manufacturing, Test and calibration of HTE vial with O2 and T spot was used The protocol for a fully automated workflow (Protocol - Operation of automated workflow for investigation of oxygen evolution) was utilized with a rare earth metal stir bar	
	- 16:35	everything was set up according to the Protocol - Operation of automated workflow for investigation of oxygen evolution	
		The python script is loaded as described in Protocol - Operation of automated workflow for investigation of oxygen evolution and the "experiment.yml" is initialized with four experiments (see reagents table for volumes of reactants) with the "run"-parameter = "true" and one at the end of the queue with "run" = "false"	
		The python script is executed with "python run.py"	
		AutoSuite program was started -->priming of tubings	
	16:40	from this point on the execution was done fully automatic (except the changing of lids and refill of stock vials) according to the table with steps executed by the automated platform in the Protocol - Operation of automated workflow for investigation of oxygen evolution	
	ca. 16:50	MRG-059-ZP-1-1 was started	septum was pierced 0 times at this point
	ca. 18:10	MRG-059-ZP-2-1 was started	at this point the setpum in the vial lid was pierced 9 times
	ca. 19:30	MRG-059-ZP-1-2 was started	at this point the setpum in the vial lid was pierced 18 times
	ca. 20:50	MRG-059-ZP-2-2 was started	at this point the setpum in the vial lid was pierced 27 times
	ca. 22:10	MRG-059-ZP-1-3 was started	at this point the setpum in the vial lid was pierced 36 times

	approx. 23:30	ZP-1-3 end	at this point the setpum in the vial lid was pierced 45 times
--	------------------	------------	---

Results

Analysis

done with hte_streamlit (<https://github.com/jschneidewind/hte-streamlit>)

Standard conditions look good.



pH = 9.2



-ZP-1-1 looks like outlier, other two are at 2.1. In comparison to -ZO-7: ok fit (-7 has rate of 2.24 and 2.33) -> should be in margin of error. -20 much higher rate (approx. 2.9) does not fit to other experiments, therefor likely double fluke or error in run.

Highlights importance of addition of known conditions in each run.

Linked experiments

- AE-262: Irradiation of PhPDA (AE-257), 1.5 mg/mg SDS, 2 mg/mL PhPDA after 3 d

- AE-265: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$
- AE-266: Calibration of oxygen and temperature sensor spots in the 10 mL HTE vial
- AE-267: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$
- AE-271: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, with manual light control
- AE-272: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, with manual light control I
- MRG-059-A: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals
- MRG-059-B: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals II
- MRG-059-C: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals III
- MRG-059-D: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals IV
- MRG-059-E: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals V
- MRG-059-F: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals VI
- MRG-059-Q: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try I
- MRG-059-R: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try II
- MRG-059-S: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try III
- MRG-059-T: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try IV
- MRG-059-U: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, using a rare earth metal stir bar, fully automated peripherals - 4th try V
- AE-JSC-321: Manufacturing and calibration of new 10 mL HTE with sensor spots I
- AE-323: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$
- MRG-059-V: Irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{ H}_2\text{O}$ in the ChemSpeed robot, fully automated peripherals

- 1st screening of Ru-cat concentration

- MRG-059-W: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O screenings (MRG-059-X, -Y, -Z)

- MRG-059-X: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated peripherals
- 1st screening of sacrificial oxidant

- MRG-059-Y: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated peripherals -
2nd screening of [Ru]

- MRG-059-Z: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated peripherals
- 1st screening of pH (9.2 - 10.4)

- MRG-059-ZA - Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - 2nd pH screening (5.8
-12)

- AE-341: Power measurment with different power settings from Joy-It with 470 nm

- AE-342: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O I

- MRG-059-ZB: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated
peripherals - screening of light intensity

HTE - MRG-059-ZI: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated
peripherals - screening of light intensity I

HTE - MRG-059-ZJ: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated
peripherals - reproduction of standard conditions - new vial (with ruthenium stock solution or directly
weighed in)

HTE - AE-379: Test and calibration of O₂ vial for O₂ and T measurment in irrad setup

HTE - AE-383: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O

HTE - MRG-AE-059-ZK: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot, fully automated
peripherals - reproduction of standard conditions - new vial, new T calibration

HTE - AE-MRG-059-ZL: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - reproduction of
standard conditions - new calibration

HTE - AE-404: Manufacturing, Test and calibration of HTE vial with O₂ and T spot

HTE - AE-405: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O I

HTE - AE-406: MRG-059-ZM: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - reproduction
of standard conditions - new vial

HTE - AE-413: Titration of Na₂CO₃ with NaHCO₃ and vice versa

HTE - AE-414: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O

HTE - MRG-059-ZN: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - screening of pH, [Ru], [Ox] and Irradiance

HTE - AE-422: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O

HTE - MRG-059-ZO: Irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O in the ChemSpeed robot - repetition of experiments from -ZN with not reproducible results

HTE - AE-434: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O

Linked items

Equipment - [Irradiation Set-Up](#)

Equipment - [Joy-it JT-RD6006 DC POWER SUPPLY](#)

Light Source - [UHP LED 470 nm](#)

Protocol - [Operation of automated workflow for investigation of oxygen evolution - as of 03.07.2024](#)

Attached files

newplot26.png

sha256: 214b00a9a03b5c3b43c1c09ba8719e8c6946ec62847728ce0a1c8985e276dcd6



newplot25.png

sha256: 36983452c79904491516ce428a658a4af036ae8e75b1b15a6e256e7f4dc371fd



MRG-059-ZP.zip

sha256: c98ff57741ba1c710b731e9085a187109fa615630e48dc0f48cc6d0911f1fd9a

AE-434-HTE.xlsx

sha256: 479edddc401588aeae7aa62c4f774069efb8adca818817665c66e8e8a9bf441f



Unique eLabID: 20250203-97248a81851ea512ea8f3ba76348b4531468afbc
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1690>