

AE-422: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O

Date: 2025-01-21


Tags: [Ru(bpy)₃]Cl₂*6 H₂O Stocksolution AE HTE

Category: HTE

Status: Done

Created by: Alexander Eith

Reaction scheme/sample structure



[Ru(bpy)₃]Cl₂ * 6H₂O
Chemical Formula: C₃₀H₃₆Cl₂N₆O₆Ru
Molecular Weight: 748,62400

Chemical Formula: CHNaO₃
Molecular Weight: 84,00577

Na₂S₂O₈
Chemical Formula: Na₂O₈S₂
Molecular Weight: 238,09154

NaHCO₃

Na₂CO₃
Chemical Formula: CNa₂O₃
Molecular Weight: 105,98754

H₂O
water
Chemical Formula: H₂O
Molecular Weight: 18,01500

Literature/reference experiments

Literature	https://doi.org/10.1021/acscatal.6b02595
Reproduction	/
Related experiment	HTE - AE-383: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O HTE - AE-405: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O I HTE - AE-414: Preparation of stock solutions for the irradiation of [Ru(bpy)₃]Cl₂ * 6 H₂O

Prepared solutions

Name	Chemical	Cas Number	Amount [mmol]	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume(added Milli-Q water) [mL]	Optained conc.
Ru II	[Ru(bpy) ₃]Cl ₂ *6 H ₂ O	50525-27-4	3.0 μmol	2.246	2.15	748.62	14.36	200 μM
Ox III	Sodium persulfate	7775-27-1	1.8	428.56	428.57	238.09	15	120 mM
NaHC O3 I	Sodiumhydrogen carbonate	144-55-8	7.2	604.84	605.26	84.006	10	0.72 M
NaHC O3 II	Sodiumhydrogen carbonate	144-55-8	15	1260.1	1260.11	84.006	15	1 M
NaHC O3 III	Sodiumhydrogen carbonate	144-55-8	0.5	42.00	42.18	84.006	10	0.05 M
Na2C O3 I	Sodium carbonate (anhydrous)	497-19-8	15	1589.8	1590.04	105.988	15	1 M

Na ₂ CO ₃ II	Sodium carbonate (anhydrous)	497-19-8	3	317.96	317.84	105.988	10	0.3 M
------------------------------------	------------------------------	----------	---	--------	--------	---------	----	-------

Procedure/observations

Date	Time	Step	Observations
21.01	9:30	Each chemical, was weighed in a separate 15 mL snap-on cap vial.	
	9:50	To each vial Milli-Q water (see Prepared solutions table) was added using a 2 figure scale.	
	9:50	All samples were shortly vortexed to obtain/ensure a homogeneous solution (Vortex VWR VV3, 4/6 intensity of shaking)	
	10:05	The solutions were transferred into Chem Speed Vials	

Linked experiments

- AE-265: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- AE-270: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- AE-311: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ I
- AE-323: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- AE-334: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- AE-342: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ I
- HTE - AE-365: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- HTE - AE-368: Preparation of stock solutions of $[\text{Ru}(\text{pby})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- HTE - AE-369: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- HTE - AE-372: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- HTE - AE-380: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ I
- HTE - AE-383: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- HTE - AE-395: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$
- HTE - AE-405: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$ I
- HTE - AE-414: Preparation of stock solutions for the irradiation of $[\text{Ru}(\text{bpy})_3]\text{Cl}_2 \cdot 6 \text{H}_2\text{O}$



Unique eLabID: 20250121-84da50531a601537dcdeea33fed2c4377094fb58
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1653>