

AE-412: Calibration of pH-Meter

Date: 2025-01-08

Tags: pH AE Calibration pH meter

Category: Organisational

Status: Done

Created by: Alexander Eith

Literature/reference experiments

Literature	/
Reproduction	/
Related experiment	/

Reagents

Name	Order number (Carl Roth)	Inventory number	Volume [ml]
pH-Pufferlösung ROTILABO® pH 4,01 in Beuteln	KAY0.1	/	25
pH-Pufferlösung ROTILABO® pH 7,00 in Beuteln	KAY1.1	/	25
pH-Pufferlösung ROTILABO® pH 10,01 in Beuteln	KAY2.1	/	25

Procedure/observations

Date	Time	Step	Observations
07.01	14:15	The 3 point calibration was done according to the manual by clicking the calc button and waiting till signals is stable for pH = 7.00, 4.00, 10.01 in that order.	box for pH meter.jpg pH meter.jpg electorde.jpg setup.jpg

Results

Displayed 100 %

Linked experiments

- [AE-344: Calibration of pH-Meter](#)

Organisational - [AE-394: Calibration of pH-Meter I](#)

Attached files

electorde.jpg

sha256: a570f6baf394ce15866bc2eabc02f545d13f608cfb9bc9f72c14ca35a6628d7d



box-for-pH-meter.jpg

sha256: 04147014d9ca22c50c17726ef8d71693b86b66315ce12e7ca9f201af2bd4fb5e



pH-meter.jpg

sha256: 0e1d53d545dd7775203913381f1435eb6eb0cd9bb42323e777ef24d709fe23af



setup.jpg

sha256: 05a326b93c4eb88408327ea2c68ba80f12dcb750a5d56a3d64b9fa1fd6e07186



Unique eLabID: 20250108-2a07ce9dcb6e9bf8d149d56fc7b9344bf2449bbe
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1612>